

KombiKon™

The world class needle free sampling connector with unique observation window for fast and easy sample taking



Safer • Easier • Faster

KombiKon™

Needle stick and infection protection

For more than 20 years our renowned KombiKon™ needle-free sampling connector has been protecting caregivers from needle stick injuries, one of the most common dangers to caregivers in their working environment¹. With its unique design it furthermore helps to reduce the contamination risk to the drainage system and thus the risk of retrograde infection².

Innovative design enhancements improve the KombiKon™ further, making it even safer, easier and more convenient to use.

Safer


Christmas Tree design with an increased number of ridges, improves catheter fit and retention

Easier

Integral retention ring, designed to hold the protective cap securely in place, making operation more convenient

Faster

Transparent window on the front side allows quick and easy observation of urine build up when clamping the collection system for sample taking



KombiKon™ is available for the following product ranges:



UnoMeter™ 500



Unometer™ Abdo-Pressure™

¹ 'Needlestick Prevention for Worker Safety: Sharps Injuries Still a Leading Occupational Hazard', Environment of Care News, August 2005, Volume 8, Issue 8.
² Performance of KombiKon needle free sample port. July 2004 Data on file, Unomedical a/s, a ConvaTec company