



# ORIGO™

NEEDLE GUIDE SYSTEM

STABILITY + FUNCTIONALITY + RELIABILITY = ORIGO™  
COST EFFICIENCY



ORIGO™ - From the maker of the original Swedish  
bolt-action needle guide design



Retrofits all Ultra-Pro™ and Ultra-Pro II™ needle guide adaptors\*

# ORIGO™

## DESCRIPTION



Unique Clip and Needle Pin Set designed with a user friendly and stable Quick-Needle-Release function

### ORIGO™ Needle Guide features and benefits

- Mount design that installs with one easy click for optimal patient safety
- The repeatable Quick Needle Release bolt-action ensures that the system is stable and will not fall apart during the procedure, and enables easy needle detachment
- Range of 10 gauge sizes 8,5F, 14-23G, makes it suitable for all kind of biopsy and drainage procedures
- Easy-to-read gauge sizes for convenient use in darkened ultrasound suites
- Minimized risk for needle deflection during procedures for optimal accuracy and performance
- The ORIGO™ needle guide kit includes Ascendia's patent pending UltraFit® latex free transducer cover ensuring a snug fit, which minimizes air bubbles and weld ears
- Compatible with all Ultra-Pro™ and Ultra-Pro II™ needle guide adaptors\*, allowing Ascendia MedTech AB to be your single-source general purpose needle guide provider

ORIGO™ Needle Guide System kits	
Article #	Description
160-1024	ORIGO™ needle guide (Guide only)
210-1624	ORIGO™ needle guide with 14 x 91,5 cm UltraFit® transducer cover
210-1824	ORIGO™ needle guide with 14 x 147 cm UltraFit® transducer cover
<b>Kit Contents:</b>	
24 procedure kits per box, each kit wrapped in surgical drape, containing:	
1Ea ORIGO™ Needle guide: needle Clip + Needle Pin Set 8,5F, 14-23G (not 19G)	
1Ea UltraFit® transducer cover with 1Ea SafeScan™ gel, 2Ea non-latex rubber bands and 2Ea fasteners	

Distributed by:



ASCENDIA MEDTECH AB  
 Finlandsgatan 18 • SE-16474 Kista • SWEDEN  
 Phone + 46-8-444 5400 • Fax +46-8-623 0972  
 Email: info@amedic.se  
 www.amedic.se